



Infrared Thermometer

[Overview]

This Infrared Thermometer, type AET-R161, is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of three months and above of people by detecting infrared heat from the center of the forehead.

[Features]








Automatic shutdown




Fever warning


[Unit Description]





[Display information]






Changing Temperature Scales


- 1.Ensure the thermometer is switched off.
- 2.Press the " \[\]" button for 8 seconds to bring it into temperature mode. The current temperature unit will be displayed.
- 3.Release and press the" \[\]" button again to toggle and select between °C and °F.
- 4.The thermometer will switch off automatically if it idles for more than 4 seconds and will display the last reading.

[Measurement method]





2. Place the thermometer probe point near the forehead and press" \[\]" and release after ~1 second and the thermometer will beep and display the accurate temperature. Whenever consecutive readings are required, please wait for the --- °C to flash. The above procedure can be repeated as many times as required.





3. To turn off the device press and hold the " \blacksquare " switch for \sim 5 seconds and the thermometer will turn off automatically after \sim 60 seconds.



[General Specifications]

Parameter	Specifications
Measurement Mode	Forehead mode
Measurement Units	Celsius(°C)and Fahrenheit(°F)
Operating Conditions	16°C-35°C(60.8°F to 95 °F) with a relative humidity of 15%-80%
Storage Conditions	-20°C-+55°C(-4°F to 131 °F) with a relative humidity of 15%-93%
Measurement distance	15mm to 50mm from the measuring point
Forehead temperature measurement range and accuracy	Measurement range: 32.0°C-42.2°C (89.6°F-108°F) Measurement accuracy: ±0.2°C (35.0°C-42.0°C)/ ±0.4°F(95.0°F-107.6°F)
Memories recall	32 measurements recall
Dimensions	150mmx 37mmx 16.3mm
Weight	51 g without batteries
Battery	2XAA batteries (DC.3V)
Switch off	Automatically switches off after 1 minute





Accessibility Contact FDA Careers FDA Basics FOIA No FEAR Act Nondiscrimination Website Policies



Acknowledgment Letter

2/5/2018

Reanny Wang, Medical Device Consultant Shenzhen Reanny Medical Devices Management Consulting Co, Lt Room 2012 of Gebu commercial building, Hongxing, Songgang Baoan District Shenzhen, Guangdong 518000 CHINA

Dear Reanny Wang:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has received your submission. This submission has been assigned the unique document control number below. All future correspondence regarding this submission should be identified prominently with the number assigned and should be submitted to the Document Control Center at the above letterhead address. Failure to do so may result in processing delays. If you believe the information identified below is incorrect, please notify the Program Operations Staff at (301) 796-5640.

Submission Number: K180207

Received: 1/24/2018

Applicant: Alicn Medical (Shenzhen), Inc Device: Non-contact Infrared Thermometer

We will notify you when the review of this submission has been completed or if any additional information is required. For information about CDRH review regulations and policies, please refer to http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm.

Sincerely yours,

Center for Devices and Radiological Health

U.S. Food & Drug Administration 10903 New Hampshire Avenue Silver Spring, MD 20993 www.fda.gov





CERTIFICATE

No. Q5 17 12 78453 013

Holder of Certificate: Alicn Medical Shenzhen, Inc

4/F, B Building

Shenfubao Modern Optical Factory, Kengzi Street

Pingshan District 518122 Shenzhen City

PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate:



Design and Development, Production and Sales of Clinical Electrical Thermometer, Clinical Infrared Thermometer, Blood Pressure Monitor and Video

Laryngoscope

The Certification Body of TÜV SÜL) Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.:

GZ1705201

Valid from: Valid until: 2018-06-04 2021-02-07

Date, 2018-06-04

Stefan Preiß

Page 1 of 2









CERTIFICATE No. Q5 17 12 78453 013

Applied Standard(s): EN ISO 13485:2016

Medical devices - Quality management systems -

Requirements for regulatory purposes

(ISO 13485:2016) DIN EN ISO 13485:2016

Facility(ies):

Alicn Medical Shenzhen, Inc.

4/F, B Building, Shenfubao Modern Optical Factory, Kengzi Street, Pingshan District, 518122 Shenzhen City, PEOPLE'S

REPUBLIC OF CHINA





