


# CE Protective Mask CE




## [Standard]

GB / T20944.3- 2008


### [Features]

- **Strong antibacterial layer**  
(Copper ions in contact with bacteria can destroy the structure of bacteria and inhibit bacterial reproduction and proliferation)
- **Antibacterial, dust-proof, anti-pollen allergy**
- **Anti-vehicle exhaust, anti-haze anti-odor, anti-PM2.5**
- **Accordance with the special technology of human facial engineering**
- **Efficient combination of fit and permeability.**

### [Material ]



# CE Protective Mask CE



[Product]: Daily Disposable Protective Mask

[Model]: 22x13.5

[Executive Standard]: GB 20944.3-2008

[Expiry Date]: 3 Years

[Scope of Application]: Using for ordinary people except infants and children to block and filter the non-oil particles in the air in normal life.

[Instructions]:



1. Take the folded mask out of the package. Gently Open the mask while nose dip pointing upwards.



2. Holding each side of the ear bands to fix the mask on the face to cover the mouth and nose.



3. Gently press the nose dip with your hand to fit the height and shape of the nose bridge. Check the fitness between face and the edge of the mask.



4. Adjust the each ear band until feeling comfortable. When the mask wearing properly, there should be a slight indentation while breathing.

[Matters For Attention]:

1. If feel uncomfortable during wearing, should stop using immediately
2. This product is not suitable for using in extreme weather, anoxic environment, underwater operation, escaping, fire fighting, special industries or other certain professions.
3. Discontinue use of the product when odor, contamination, dampness, or discoloration occur.
4. People allergic to nonwoven fabrics DO NOT use.
5. The product should be inspected for damage and contamination before use.
6. This product is not designed for wash.



# EU Declaration of Conformity

**Manufacturer:** Chengdu Medical and Sanitary Material Factory  
Chengdu Economic and Technological Development Zone, Sichuan  
Province (14km of The old Chengdu-Chongqing Highway), The People's  
Republic of China.

**European Representative:** Luxus Lebenswelt GmbH  
Kochstr.1, 47877, Willich, Germany

**Product Name:** Single-use Medical Face Mask  
**Specification:** Flat/3D  
**GMDN Code:** 35177  
**UMDN Code:** 12-458

**Classification (MDR, Annex VIII):** Class I, Rule 1.

**Conformity Assessment Route:** EU DECLARATION OF CONFORMITY following the Annex II +  
Annex III + Article 19 of MDR (EU) 2017/745.

We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EU Regulation and Standards. All supporting documentations are retained under the premises of the manufacturer.

Chengdu Medical and Sanitary Material Factory is exclusively responsible for the declaration of conformity.

**General applicable regulations, directives:**

Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC.

**Applied standards, common specification, guidance:**

EN 14683:2019+AC:2019, EN ISO 15223-1:2016, EN 1041:2008, EN ISO 14971:2012, EN 62366-1:2015+AC:2015, ISO 10993-1:2018, ISO 10993-5:2009, ISO 10993-10:2010, ASTM D4169-2016, MDCG 2019-15.

SIGNATURE

Name: Zeng Qunfang  
Position: Production Manager  
Place/date: 25th. Mar. 2020



CE